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
Joseph Klebba, Ph.D

"I have two young children; I am no longer familiar with the concept of free time." Nevertheless, Simple Western's hands-free automation and rapid time to results have given him the flexibility to spend his time working on other experiments or spending time with his family.


From your peers: Targeting Once Inaccessible Phospho-Proteins with RePlex :: ProteinSimple

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実例：Professor Dr. Rieger

**LOEWE Center
for Cell & Gene
Therapy**
Goethe University
Hospital Frankfurt,
Germany

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
The screenshot shows the 'New Jess Assay' configuration window. On the left, the 'RePlex' assay is selected, and 'Chemiluminescence' is chosen for both Probe 1 and Probe 2. The '2 IMMUNOASSAYS' plate layout is shown on the right, with wells A-K and 1-25. Red boxes highlight the assay selection, the detection method, and the reagent allocation for the second assay (E, F, J, K rows).

Legend for Plate Layout:


- Row A:** Biotinylated Ladder, 5 μL ; Prepared Samples, 3 μL
- Row B:** Antibody Diluent (2 or Milk-Free), 10 μL
- Row C:** Antibody Diluent (2 or Milk-Free), 10 μL ; Primary Antibody for Probe 1, 10 μL
- Row D:** Streptavidin-HRP or NIR, 10 μL ; Secondary Antibody for Probe 1, 10 μL
- Row E:** Antibody Diluent (2 or Milk-Free), 10 μL ; Primary Antibody for Probe 2, 10 μL
- Row F:** Antibody Diluent (2 or Milk-Free), 10 μL ; Secondary Antibody for Probe 2, 10 μL
- Row G:** Wash Buffer
- Row H:** Wash Buffer 500 μL /compartment
- Row I:** Wash Buffer
- Row J:** Luminol-Peroxide Mix, 170 μL /compartment (for chemiluminescence only)
- Row K:** RePlex[®] reagent mix, 300 μL /compartment

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


Total Protein

Reagent List:


- Biotinylated Ladder, 5 μL ; Prepared Samples, 3 μL
- Antibody Diluent (2 or Milk-Free), 10 μL ; Total Protein labeling reagent, 10 μL
- Antibody Diluent (2 or Milk-Free), 10 μL
- Antibody Diluent (2 or Milk-Free), 10 μL ; Primary Antibody for Probe 1, 10 μL
- Streptavidin-HRP or NIR, 10 μL ; Secondary Antibody for Probe 1, 10 μL
- Total Protein Streptavidin-HRP for Probe 2, 8 μL
- Wash Buffer
- 500 μL /compartment
- Luminol-Peroxide Mix, 170 μL /compartment
- RePlex™ reagent mix, 300 μL /compartment

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
Case 3 : Housekeeping proteins *versus* Total Protein

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«Cell confluence significantly affects the levels of α -tubulin and Glyceraldehyde-3-Phosphate Dehydrogenase (GAPDH)

Levels of heat-shock protein-90 (α and β) and β -actin remained unchanged at a wide range of cell densities»

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
[Housekeeping genes; expression levels may change with density of cultured cells - ScienceDirect](#)

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
«Housekeeping» protein signal depends on cell density

RePlex Webinar:

[The Benefits of Sequential, Multiplexed Westerns using RePlex on Simple Western: 3 case studies from University of Lausanne | ProteinSimple](#)

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
Fig. S1. Western blotting for dystrophin using capillary immunoelectrophoresis with total protein assay. A) WT mouse (cranialis tibialis) and dog (vastus lateralis) skeletal muscle express only full-length dystrophin, confirming that the band detected at ~280kDa is Dp427, as has been previously reported and likely a consequence of the molecular weight ladder used, which underestimates molecular weights above 280 kDa (Beekman et al., 2018).

B) Total protein loading (TP) per lane is shown adjacent to dystrophin labelling of each sample. C) Dystrophin protein expression was evaluated via capillary immunoelectrophoresis (ProteinSimple Jess), using a C-terminal antibody (AbCam154168, 1:250) and anti-mouse HRP (ProteinSimple). Lysates were loaded at $0.1\mu\text{g}\cdot\mu\text{l}^{-1}$, and all dystrophin signal intensities were normalized to total protein load (determined via Total Protein RePlex assay as per the manufacturer's instructions). Relative abundance was calculated for the full data set by using the maximum value across all isoforms.

Dystrophin protein expression was evaluated via capillary immunoelectrophoresis (ProteinSimple Jess), using a C-terminal antibody (AbCam154168, 1:250) and anti-mouse HRP (ProteinSimple). Lysates were loaded at $0.1\mu\text{g}\cdot\mu\text{l}^{-1}$, and all dystrophin signal intensities were normalized to total protein load (determined via Total Protein RePlex assay as per the manufacturer's instructions). Relative abundance was calculated for the full data set by using the maximum value across all isoforms.

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OK

● Biotinylated Ladder, 5 μ L ● 希釈済みサンプル, 3 μ L

● Milk-Free Antibody Diluent, 10 μ L ● Protein Normalization Reagent, 10 μ L

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Protein Normalization Module (DM-PN02)使用

a biotechnne brand


Figure S7. Immunoblot analysis of the expression level of three liver RP proteins show no changes. (A) Representative immunoblot of polysomal fraction from liver samples analyzed by capillary electrophoresis and probed with antibodies recognizing eL22 (o), RACK1 (*) and uS3 (<). The immunofluorescence signal is overlaid with signal from the total protein staining. Protein weight markers in kDa are shown on the left of each blot. (B) The signal for each RP from the immunoblot (IB) is normalized to the total protein concentration of each samples using the total protein staining. Data are means \pm SEM, n=4-6 biological replicates. For comparison, the RP expression level of these three proteins from the LC/MS-MS are

Immunoblot Analysis. The immunoblot analysis was performed with liver samples using the capillary electrophoresis system (Jess, ProteinSimple). Frozen tissues (n = 4 to 6 biological replicates) were homogenized separately and fractionated in monosomal and polysomal fractions as prepared for the mass spectrometry analysis (see above) and loaded on the Jess separation module (2 to 40 kDa and 12 to 230 kDa). Antibodies against RACK1, uS3, and eL22 were purchased from Santa Cruz Biotechnology. The corresponding RP peaks were normalized to the total protein concentration using a protein normalization kit (ProteinSimple) and quantified using the Jess quantification module.

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[Invariable stoichiometry of ribosomal proteins in mouse brain tissues with aging | PNAS](#)

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


REAGENTS	NO MULTIPLEX (ONE COLOR)	MULTIPLEX (TWO COLORS)
Stellar Total Protein Streptavidin-HRP	288 µL	276 µL
Secondary Antibody	6 µL (Rabbit or Mouse)	6 µL Rabbit 6 µL Mouse
Oligo	6 µL (Rabbit or Mouse)	6 µL Rabbit 6 µL Mouse
Total Volume	300 µL	300 µL

- Vortex mixture for 30 seconds to mix thoroughly.
- Combine 450 µL Luminol-S and 450 µL Peroxide in a microcentrifuge tube.

- The Biotin Labeling Reagent should only be prepared immediately before loading the plate.

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
Figure 10. Western blot analysis of mitochondria and ER-associated proteins. (a) Whole cell lysates and mitochondrial lysates of parental CHO and MMP-enriched host cells were analyzed for proteins involved in mitochondrial function and high MMP phenotype. (b and c) Whole cell lysates of stable pools generated from parental CHO and MMP-enriched hosts were analyzed for mitochondrial and ER-stress proteins. Beta-actin and succinate dehydrogenase subunit A were used as loading controls. Western blot analysis of the lysates of host cells and pools demonstrating upregulation of proteins associated with mitochondrial functions and downregulation of ER-stress related proteins in MMP-enriched host cells and pools.

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[Mitochondrial membrane potential-enriched CHO host: a novel and powerful tool for improving biomanufacturing capability \(nih.gov\)](https://doi.org/10.1038/nbt.3588)

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
Fig. 7. Induction of apoptotic proteins caused by pesticides mixtures at concentration of 30 μ M. a: The lane graph of apoptotic protein mixture, DC stands for difenoconazole + cypermethrin, DT stands for difenoconazole + triazophos, Con stands for the control. b: The apoptotic protein content in HepG2 cells compared with the control.

Note: * and ** indicate that the difference is significant ($P < 0.05$) and highly significant ($P < 0.01$) compared with the difenoc

2.8. Western blot analysis using the Wes-ProteinSimple system

Based on previous results, some pesticide mixtures with obviously reduction of cell viability were selected to evaluate protein levels in HepG2 cells after 24 h exposure at concentration of 30 μ M. Then the treated cells were lysed on ice with the cell lysis buffer for Western and immunol precipitation (IP) containing protease and phosphatase inhibitor (Chen et al., 2010; Zhang et al., 2007). The protein content was measured by bicinchoninic acid (BCA) protein assay kit (Walker, 1994). Quantitative protein analyses were performed according to the ProteinSimple user manual (Chen et al., 2013; Yu et al., 2020). In brief, cell lysate samples were mixed with the loading buffer and dithiothreitol (DTT), and then heated at 95 $^{\circ}$ C for 5 min. The samples, fluorescent molecular weight markers, primary antibodies, HRP-conjugated secondary antibodies, chemiluminescent substrate, and separation and stacking matrices were dispensed in a designated plate. Following plate loading, the separation electrophoresis and immunodetection was conducted automatically in the capillary system. The digital image was analyzed with Compass software (ProteinSimple, USA), and the quantified data of the detected protein were reported as molecular weight, signal/peak intensity. At least three independent experiments were conducted and the results were expressed as the percentage of protein content (%) with respect to the control.

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[Three widely used pesticides and their mixtures induced cytotoxicity and apoptosis through the ROS-related caspase pathway in HepG2 cells - ScienceDirect](#)

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