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
Joseph Klebba, Ph.D

"I have two young children; I am no longer familiar with the concept of free time." Nevertheless, Simple Western's hands-free automation and rapid time to results have given him the flexibility to spend his time working on other experiments or spending time with his family.

From your peers: Targeting Once Inaccessible Phospho-Proteins with RePlex :: ProteinSimple

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
[Human adipose tissue protein analyses using capillary western blot technology | Nutrition & Diabetes \(nature.com\)](#)

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
[Development of an Automated Capillary Immunoassay to Detect Prion Glycotypes in Creutzfeldt-Jakob Disease - Laboratory Investigation \(uscap.org\)](http://uscap.org)

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

FIGURE 3. Simple Western analysis of HQC, MQC, and LQC samples.

HQC, MQC, and LQC samples contain rheNOS at final concentrations of 15, 3, and 0.6 ng/mL, respectively. Each sample was probed with the anti-eNOS antibody. Each QC sample (n=2) was subtracted by the matrix-only control and then the percent recovery of each sample was calculated and plotted against expected rheNOS concentrations.

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実例 : Professor Dr. Rieger

LOEWE Center
for Cell & Gene
Therapy
Goethe University
Hospital Frankfurt,
Germany

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




protein simple[®]

a biotechne brand

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Simple Western Antibody Database | Bio-Techne

Simple Western Validated Antibodies

Filter Antibodies: (3864)

| Protein Target | Antibody Type | Vendor / Cat # | Protein Isoform | Citations & Supporting Data | Host Species | Sample | Separation | Dilution | kDa on Simple Western |
|----------------|---------------|----------------------------------|---------------------------|-----------------------------|-------------------|--|------------|----------|-----------------------|
| α-Nucleocapsid | Primary | Novus Biologicals NB100-56683 | SARS Nucleocapsid Protein | Datasheet | Rabbit Polyclonal | Recombinant SARS-CoV-2 Nucleocapsid Protein (Novus Biologicals PN 10474-CV); Lysates of SARS-CoV-2 | Size | 25 mg/mL | 60 kDa |
| 14-3-3 γ | Primary | Novus Biologicals NB100-406 | 14-3-3 gamma (γ) | Datasheet | Mouse Monoclonal | HeLa lysate 1.0 mg/ml | Size | 1:12.5 | 36 kDa |

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Instrument Citation Search Tool | Bio-Techne

Instrument Citation Search

At Bio-Techne, we're changing the way scientists analyze proteins. Our innovative product portfolio helps researchers reveal new insight into proteins, advancing their understanding of protein function. As a result, our products are rapidly increasing in the scientific literature across many applications. Use this citation database to find scientific articles that feature our instruments.

Showing 1 - 25 of 1407 citations

Filter by:

Instrument Platform Research Areas Year Clear filters X

Search Title:

Instrument Platform:

Simple Plex

Research Category:

Cell Therapy

Bourbon E et al. (2022) HLA-DR expression on monocytes and outcome of anti-CD19 CAR-T cell therapy for large B-cell lymphoma Blood Adv.

PMID: 35439292

Instrument Platform:


Simple Plex

Research Category:

COVID-19

Fabris M et al. (2022) High T-cell response rate after COVID-19 vaccination in belimumab and rituximab recipients J Autoimmun.

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
The screenshot shows the 'New Jess Assay' configuration window. On the left, the 'RePlex' assay is selected, and 'Chemiluminescence' is chosen for both Probe 1 and Probe 2. The '2 IMMUNOASSAYS' plate layout is shown on the right, with wells A-K and 1-25. A legend on the right lists the reagents for each well, with red boxes highlighting the reagents for E and F rows, and the reagents for J and K rows.


Legend:

- Biotinylated Ladder, 5 μL ; ● Prepared Samples, 3 μL
- Antibody Diluent (2 or Milk-Free), 10 μL
- Antibody Diluent (2 or Milk-Free), 10 μL ; ● Primary Antibody for Probe 1, 10 μL
- Streptavidin-HRP or NIR, 10 μL ; ● Secondary Antibody for Probe 1, 10 μL
- Antibody Diluent (2 or Milk-Free), 10 μL ; ● Primary Antibody for Probe 2, 10 μL
- Antibody Diluent (2 or Milk-Free), 10 μL ; ● Secondary Antibody for Probe 2, 10 μL
- Wash Buffer 500 μL /compartment
- Luminol-Peroxide Mix, 170 μL /compartment (for chemiluminescence only)
- RePlex[®] reagent mix, 300 μL /compartment

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



2.8. Simple Western analysis

For analysis of PKA, total-ERK1/2, and ratio of p-ERK1/2/Total ERK1/2 protein abundance, simple Western analyses were performed using the RePlex method for Jess capillary Western blot system (Protein Simple, San Jose, Ca) following the manufacturer's protocol. The total-ERK1/ERK2 and p-ERK1/ERK2 relative abundance measurement were performed in the same plate, but in separated wells using 0.108 µg per reaction. PKA relative abundance analyses were performed using 0.216 µg per reaction. Briefly, samples were mixed with a master mix (Protein Simple, San Jose, Ca), heated at 95 °C for 5 min, and 3 µL of each sample were loaded. Each capillary from the 12–230 kDa cartridge (Protein Simple, cat#SM-W004) aspirated the separation matrix followed by the stacking matrix and the protein lysates. Capillaries were lowered to be in contact to running buffer and voltage was applied for appropriate protein separation by molecular weight. Following, UV light immobilized the protein to the capillary wall and target proteins were immunoprobed with primary antibodies (ERK1/2, Novus Biological cat # NBP2-12728; p-ERK1/ERK2, R&D Systems cat#AF1018; or PKAc-α, Cell Signaling, cat# 4782S) diluted in antibody diluent 2 buffer (Protein Simple, cat#042-203) followed by HRlabeled secondary antibody for improved chemiluminescent detection. For each analysis, the target protein antibody from the first immunoassay was removed using the RePlex reagent (Protein Simple, cat#RP-001) and a second immunoassay was performed in the same capillary using the anti-BActin antibody (R&D Systems, cat#MAB3929) and used for standardization. Data were generated using Compass for SW 6.0 for Windows (Protein Simple, San Jose, Ca). For data analysis, the area of spectra that matched the target protein was used.

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[Effects of asprosin on estradiol and progesterone secretion and proliferation of bovine granulosa cells - ScienceDirect](#)

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

Case 3 : Housekeeping proteins *versus* Total Protein

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«Cell confluence significantly affects the levels of α -tubulin and Glyceraldehyde-3-Phosphate Dehydrogenase (GAPDH)

Levels of **heat-shock protein-90** (α and β) and **β -actin** remained unchanged at a wide range of cell densities»

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
[Housekeeping genes; expression levels may change with density of cultured cells - ScienceDirect](#)

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
«Housekeeping» protein signal depends on cell density

RePlex Webinar:

[The Benefits of Sequential, Multiplexed Westerns using RePlex on Simple Western: 3 case studies from University of Lausanne | ProteinSimple](#)

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
Fig. S1. Western blotting for dystrophin using capillary immunoelectrophoresis with total protein assay. A) WT mouse (cranialis tibialis) and dog (vastus lateralis) skeletal muscle express only full-length dystrophin, confirming that the band detected at ~280kDa is Dp427, as has been previously reported and likely a consequence of the molecular weight ladder used, which underestimates molecular weights above 280 kDa (Beekman et al., 2018).

B) Total protein loading (TP) per lane is shown adjacent to dystrophin labelling of each sample. C) Dystrophin protein expression was evaluated via capillary immunoelectrophoresis (ProteinSimple Jess), using a C-terminal antibody (AbCam154168, 1:250) and anti-mouse HRP (ProteinSimple). Lysates were loaded at $0.1\mu\text{g}\cdot\mu\text{l}^{-1}$, and all dystrophin signal intensities were normalized to total protein load (determined via Total Protein RePlex assay as per the manufacturer's instructions). Relative abundance was calculated for the full data set by using the maximum value across all isoforms.

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OK

● Biotinylated Ladder, 5 μ L ● 希釈済みサンプル, 3 μ L

● Milk-Free Antibody Diluent, 10 μ L ● Protein Normalization Reagent, 10 μ L

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
a biotechnie brand

Supplemental Figure 4: Validation of Proximity Ligation Assay (PLA) technique in endothelial cells with stable ER α expression. (A) Representative ER α (green) and DAPI (blue) staining on control endothelial cells (TeloHAECs) or endothelial cells expressing ER α (ER α -TeloHAECs); Obj: X40. (B) Representative ER α protein expression and (C) total proteins analysed by Simple Western in TeloHAECs or ER α -TeloHAECs. (D) Representative images of PLA performed with ER α and SRC antibodies or with single antibodies in ER α -TeloHAECs incubated with E2 for 5min. (E) Representative image of PLA performed with ER α and SRC antibodies in TeloHAECs incubated with E2 for 5min. The detected dimers are represented by red dots. Nuclei were counterstained with DAPI (blue) (Scale bar: 20 μ m). (F) Quantification of the number of dots per cell. The experiments were reproduced 2 times.

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[NRXN1 depletion in the medial prefrontal cortex induces anxiety-like behaviors and abnormal social phenotypes along with impaired neurite outgrowth in rat | Journal of Neurodevelopmental Disorders | Full Text \(biomedcentral.com\)](#)

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

| REAGENTS | NO MULTIPLEX (ONE COLOR) | MULTIPLEX (TWO COLORS) |
|---|--------------------------|---------------------------|
| Stellar Total Protein Streptavidin-HRP | 288 µL | 276 µL |
| Secondary Antibody | 6 µL (Rabbit or Mouse) | 6 µL Rabbit 6 µL Mouse |
| Oligo | 6 µL (Rabbit or Mouse) | 6 µL Rabbit 6 µL Mouse |
| Total Volume | 300 µL | 300 µL |

- Vortex mixture for 30 seconds to mix thoroughly.
- Combine 450 µL Luminol-S and 450 µL Peroxide in a microcentrifuge tube.

• The Biotin Labeling Reagent should only be prepared immediately before loading the plate.

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
Figure 10. Western blot analysis of mitochondria and ER-associated proteins. (a) Whole cell lysates and mitochondrial lysates of parental CHO and MMP-enriched host cells were analyzed for proteins involved in mitochondrial function and high MMP phenotype. (b and c) Whole cell lysates of stable pools generated from parental CHO and MMP-enriched hosts were analyzed for mitochondrial and ER-stress proteins. Beta-actin and succinate dehydrogenase subunit A were used as loading controls. Western blot analysis of the lysates of host cells and pools demonstrating upregulation of proteins associated with mitochondrial functions and downregulation of ER-stress related proteins in MMP-enriched host cells and pools.

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[Mitochondrial membrane potential-enriched CHO host: a novel and powerful tool for improving biomanufacturing capability \(nih.gov\)](https://doi.org/10.1038/nbt.3488)

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
Fig. 5. XPO1 protein degradation is regulated by ASB8. a-b Sensitization profile of HAP1 cells with either knockout or overexpression of ASB8 under selection with selinexor. Cell viability when transduced with a safe-targeting or ASB8 targeting sgRNA (a). Cell viability when trans-

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RB, p53, Cyclin D1, p16, HPV16E7


Figure 8. Wes protein gel for extracellular vesicles and whole cell lysates from HNSCC, normal keratinocyte, and transformed cell lines. RB detected at 112 kDa, 1:100 antibody dilution, 0.5 mg/mL protein; p53 detected at 60 kDa, 52 kDa, and 48 kDa, 1:50 antibody dilution, 0.25 µg/µL protein; CyclinD1 detected at 40 kDa, 1:200 antibody dilution, 0.25 µg/µL protein; p16 detected at 24 kDa, no antibody dilution, 0.25 µg/µL protein; and HPV16E7 detected at 18 kDa, 1:100 antibody dilution, 1 µg/µL protein. Bold italics: HPV-positive.

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[Cancers | Free Full-Text | Differences in Extracellular Vesicle Protein Cargo Are Dependent on Head and Neck Squamous Cell Carcinoma Cell of Origin and Human Papillomavirus Status \(mdpi.com\)](https://www.mdpi.com/1918-4746/11/1/1)

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