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
Normalize your data with total protein detection in the same capillary using RePlex. HeLa lysates (0.01 - 0.2 mg/mL) were separated and probed with a 14-3-3 antibody in Probe 1 followed by the Total Protein Assay in Probe 2. Protein normalization using RePlex shows comparable 14-3-3 signal across increasing concentrations of HeLa lysate.

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


JBC

JOURNAL OF
BIOLOGICAL
CHEMISTRY

Record how data were obtained, whether signal intensity was linear with antigen loading, and how protein loading was normalized. Some detection methods (e.g., ECL) have a very limited linear range.

Normalize signal intensity to **total protein loading** (assessed by staining membranes for total protein) whenever possible. "House-keeping" proteins should not be used for normalization without evidence that manipulations do not affect expression. Phospho-specific antibody signals should be normalized to total levels of the target protein.

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





Absolute quantitation of endogenous proteins with precision and accuracy using a capillary Western system

Jin-Qiu Chen^{a,*}, Madeleine R. Heldman^a, Michelle A. Herrmann^a, Noemi Kedei^b, Wonhee Woo^b, Peter M. Blumberg^b, Paul K. Goldsmith^c

^a Collaborative Protein Technology Resource, Laboratory of Cell Biology, National Cancer Institute, National Institutes of Health, Bethesda, MD 20892, USA
^b Laboratory of Cancer Biology and Genetics, National Cancer Institute, National Institutes of Health, Bethesda, MD 20892, USA
^c Office of Science and Technology Partnership, Office of the Director, Center of Cancer Research, National Cancer Institute, National Institutes of Health, Bethesda, MD 20892, USA

ARTICLE INFO

Article history:
Received 16 May 2013
Received in revised form 13 July 2013
Accepted 16 July 2013
Available online 26 July 2013

Keywords:
Protein kinase C (PKC)
Capillary immunosay
Simple Western
Quantitative protein analysis

ABSTRACT

Precise and accurate quantification of protein expression levels in a complex biological setting is challenging. Here, we describe a method for absolute quantitation of endogenous proteins in cell lysates using an automated capillary immunoassay system, the size-based Simple Western system (recently developed by ProteinSimple). The method was able to accurately measure the absolute amounts of target proteins at picogram or sub-picogram levels per nanogram of cell lysates. The measurements were independent of the cell matrix or the cell lysis buffer and were not affected by different antibody affinities for their specific epitopes. We then applied this method to quantitate absolute levels of expression of protein kinase C (PKC) isoforms in LNCaP and U937 cells, two cell lines used extensively for probing the downstream biological responses to PKC targeted ligands. Our absolute quantitation confirmed the predominance of PKC δ in both cells, supporting the important functional role of this PKC isoform in these cell lines. The method described here provides an approach to accurately quantitate levels of protein expression and correlate protein level with function. In addition to enhanced accuracy relative to conventional Western analysis, it circumvents the distortions inherent in comparison with signal intensities from different antibodies with different affinities.

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
Articles:
A global transcriptional network connecting noncoding mutations to changes in tumor gene expression

Science
MAAS

Nucleic Acids Research

PNAS

PNAS
Invariable stoichiometry of ribosomal proteins in mouse brain tissues with aging

Cell
Divergent Routes toward Wnt and R-spondin Niche Independence during Human Gastric Carcinogenesis

Scientific Reports
Evaluation of ATM heterozygous mutations underlying individual differences in radiosensitivity using genome editing in human cultured cells

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Bio-Techne antibodies are now Simple Western Certified

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Human/Mouse/Rat Phospho-ERK1(T202/Y204)/ERK2 (T185/Y187) Antibody AF1018
See All Sizes & Prices ▼
10 Citations
Clonality: Polyclonal
Host: Rabbit IgG
Applications: WB, Flow, Simple Western, IHC
compare (select up to 7 products from this page to compare)

Cell Signaling TECHNOLOGY
Search by keyword or product number

- Immunohistochemistry (Frozen) (2)
- IP (3335)
- Neutralizing Ab (2)
- Peptide ELISA (34)
- RNA Dot Blot (2)
- Simple Western™ (82)
- Western blot (7132)

2933	Bim (C34C5) Rabbit mAb	WB
8547	Btk (D3H5) Rabbit mAb	WB
3195	E-Cadherin (24E10) Rabbit mAb	WB
9664	Cleaved Caspase-3 (Asp175) (5A1E) Rabbit mAb	WB
9661	Cleaved Caspase-3 (Asp175) Antibody	WB
9662	Caspase-3 Antibody	WB
14220	Caspase-3 (D3R6Y) Rabbit mAb	WB

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