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Please make the most of the facility. The mission of a shared equipment facility like ours is to support and advance your research.

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To avoid liquid nitrogen(LN) accidents, remember these three points.

- 爆発 Explosion (LN expands about 700-fold when it vaporizes.)

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Use only dewar flask designed for LN. Do not use ordinary thrmos bottle.

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A single breath of oxygen-free gas can cause you to lose consciousness.

換気が必須。Ventilation is essential.

L階液体窒素室に入室中はドアを閉めない。Keep the door open at all times while you are in the LN.

エレベータで運ぶ時は同乗しない。When transporting LN, do not ride in the elevator with it.

容器を倒さないように運ぶ。Carry the dewar flask carefully to prevent it from tipping over.

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Wear leather gloves. Do not use cotton work gloves. LN trapped in clothing causes frostbite.

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L階液体窒素室 Liquid nitrogen room at floor L

cryopreservation tank to store the samples (cells, tissues, etc)

LN fill station for laboratory use

Always enter the room in pairs.

Turn on the ventilation switch before entering the room. This will also turn on the room lights.

Keep the door open while you are in the room

Wear leather gloves.

Do not let the corners of the sample rack strike the tank.

This can damage the tank irreversibly.

**Do not ride in the elevator with LN.
Turn the warning light on.**

- See LN page in CentMeRE web <https://www.med.osaka-u.ac.jp/pub/ctrlab/ln/>
- Demonstration of liquid nitrogen on this Friday, 13:15-13:35

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After the cooling system for valuable samples failed, LN was released into the room to reduce the room temperature. Two fatalities due to hypoxia.

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For some reason, the storage tank became sealed, and the internal pressure gradually increased until it ruptured. The explosion partially destroyed the factory, and 25 buildings within a radius of 400 meters, as well as 39 parked vehicles, were damaged.

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
An case at CentMeRE

Even a small hole in the neck tube can result in loss of vacuum insulation, leading to evaporation of the LN and a rise in the tank temperature.

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The most common incidents involve placing plastic products, such as centrifuge tubes, in a dry heat sterilizer.

事故事例 An incident case at our facility

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A user attempted to sterilize 1.5 mL polypropylene tubes and containers made of polypropylene and low-density polyethylene at 200°C. These materials had a heat tolerance of 100–140°C and an ignition point of 201°C. Another user noticed white smoke coming from inside the equipment, allowing the fire to be extinguished at an early stage.

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Dry heat sterilizer during fire extinguishing.

Flames were visible immediately after the door was opened.

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P1, P2, and P3 laboratories are available for use.

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If you conduct recombinant DNA experiments at CentMeRE, you must include this in your recombinant DNA experiment application. **Notification to CentMeRE is also required.** Please check with your supervisor.

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Please use these laboratories in compliance with all applicable laws and regulations.

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Please take any liquid waste generated while using equipment at CentMeRE back to your laboratory for proper disposal.

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Although some rooms at CentMeRE have sinks used by the staff, these sinks are not available for users.

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Please feel free to ask the technical staff if you need anything.

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Shared equipment is for everyone. Please use it considerately so that all users can work comfortably. That also helps make the facility more comfortable for you.

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

If you notice anything unusual in the operation of the equipment, please call the staff. Outside business hours, please email us at staff@ctrlab.med.osaka-u.ac.jp.

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If you damage any equipment or parts, be sure to contact CentMeRE. Students will not be asked to pay for the damage.

The most serious problems arise when equipment remains out of service for a long time, or when users continue to use it while it is still not functioning properly.

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When publishing research results obtained using equipment at CentMeRE, please acknowledge the support of CentMeRE and the staff member in charge. Please also inform us once the paper has been published.

<Example> This study was supported by ([staff name] and) Center for Medical Research and Education, Graduate School of Medicine, The University of Osaka.

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Acknowledgments help make the contributions of shared equipment facility staff to your research visible. They also contribute to the introduction and renewal of equipment. (Please also acknowledge other facilities where appropriate.)

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**The staff of the shared equipment
facility are your research partners.**

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**Let's make the most of shared
equipment and do great research!**