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
---

04-21-2026

# TODAY, IDT IS THE LEADING BRAND IN DNA AND RNA SYNTHESIS

- Founded in 1987 by Dr. Joseph Walder
- Largest custom oligonucleotide manufacturer worldwide
- >1500 employees in 9 locations
- >130,000 active customers
- >95% of ordered products are manufactured and shipped in less than 24 hours


# THE WORLD'S LARGEST SUPPLIER OF CUSTOM NUCLEIC ACIDS

>64,000 oligos synthesized every day

- Major R&D teams (~80 employees)
- CRISPR
- NGS
- Genotyping
- qPCR
- Synthetic Biology
- RNAi
- Bioinformatics

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# Alt-R™ CRISPR SYSTEM: A COMPLETE WORKFLOW

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Analyze

## Alt-R Cas9 gRNA design tool

- Predesigned guides
- Custom designs
- Design checking

## Alt-R Cas9 HDR design tool

- Friendly UI
- Empirically defined design rules
- Integration with Cas9 gRNA designs

## rhAmpSeq design tool

### Alt-R gRNAs

- Cas9 crRNA:tracrRNA
- Fluorescently labeled tracrRNAs
- Cas9 sgRNA
- Cas12a crRNA
- Custom ordering for any gRNA (e.g. pegRNA, Cas13)

### Alt-R CRISPR proteins

- WT Cas9
- HiFi Cas9
- Cas9 nickases
- A.s. Cas12a *Ultra*
- L.b. Cas12a *Ultra*
- Fluorescently labeled Cas9

### Alt-R Electroporation Enhancers

### Alt-R HDR Donor Oligos

- Up to 200 nt
- Modified ssODNs

### Alt-R HDR Donor Blocks

- Modified to reduce blunt integration
- Up to 3000 nt
- Sequence-verified via NGS

### Alt-R HDR Enhancer

- V2

### Alt-R HDR Enhancer

- Protein


## rhAmpSeq system for CRISPR

- Multiplexed amplicon sequencing


## CRISPRAltRations NGS analysis tool

- Cloud-hosted UI for analysis of CRISPR on- and off-target editing

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

# Alt-R Cas9 CRISPR LIBRARIES

## Alt-R™ CUSTOM CRISPR gRNA LIBRARIES

Chemically modified guide RNA libraries for CRISPR screening research

Fast turnaround time, flexible formulation

Adaptable for different CRISPR systems

Complete CRISPR workflow solution from custom design to analysis

Alt-R Custom CRISPR gRNA libraries are available for all CRISPR nucleases, including Cas9, Cas12a, Cas13, prime editing enzymes, and others. These libraries were developed to address the need for better CRISPR screening solutions. They are chemically modified guide RNAs (gRNAs) synthesized on the IDT proprietary high-fidelity RNA manufacturing platform to provide high quality, reliable gRNA libraries with fast delivery.

### BENEFITS

1. Innovative solution offered by a global leader in RNA synthesis and CRISPR innovation
2. Reliable, consistent, and fast delivery, with custom formulations available to suit a variety of project needs
3. Adaptable for alternative CRISPR systems such as Cas12a, Cas13, and prime editing
4. Enhanced nuclease resistance for maximal editing in Cas9-expressing cells or via ribonucleoproteins (RNP)
5. Optimized RNA synthesis processes to mitigate cross-contamination risk

### PRODUCT SPECIFICATIONS

Features	Options
Design	Pre-designed, custom, user-provided
CRISPR systems	Cas9, Cas12a, Cas13, prime editing, and other alternative systems
Guaranteed yield	0.5 nmol, 2 nmol, 5 nmol, and custom normalized deliverables
Cas9 gRNA formats	2-part crRNA:tracrRNA complex and sgRNA
Custom lengths supported	30-150 nt
Chemical modifications	2'-O-methyl RNA, PS linkages, end-blocking Alt-R modifications
Plate types	96- & 384-well PCR, Deep-well, V-bottom, ECHO, custom options available
Formulation options	Multi-guide per well, pooled by gene Arrayed (single gRNA/well) Custom formulations upon request
QC	Individual ESI/MS
Supporting reagents & functional analysis pipeline (optional)	WT Cas9, HiFi Cas9, Cas12a, and Cas12a Ultra Glycerol-free options available in tubes or plates (ideal for robotics) Electroporation Enhancers rhAmpSeq™ CRISPR Analysis System (NGS-based on-/off-target editing analysis)

> [WWW.IDTDNA.COM](http://WWW.IDTDNA.COM)


### Custom gRNA libraries:

- 2 part (crRNA+tracrRNA) and sgRNA
- Standard offering provides desalted gRNAs, but purified versions, with custom mod patterns, large scale deliverables and likewise gRNAs for alternative applications built on Cas9 (i.e. Prime Editing) can easily be accommodated
- **Scales** – 2nmol, 10nmol, 50nmol, 100nmol
- Available in tube or plates
- **Custom formats** – IDT can support any format to suit your needs
  - Arrayed – individual guides targeting different genes in different wells
  - Pooled – multi-guide pooled by gene
- **Custom Plates** – IDT has a variety for you to choose from including ECHO plates
- **Need help designing?** We can assist.

# CRISPR SCREENS WITH CUSTOM Alt-R gRNA LIBRARIES

4 μM RNP Nucleofector™ Nucleofection (Lonza) into HAP1 cells, 4 gRNAs for 5 genes as 2-part or sgRNA gave >85% avg editing at all target loci

# ALT-R™ CRISPR CUSTOM GUIDE RNAS

## CUSTOM ORDERING TOOL FOR ANY gRNAs

**ORDER MENU** You are using the IDT KK portal.

### Custom CRISPR gRNA ordering

Select All ACTIONS: ▾ # of Items: 1 GO BULK INPUT 📄

# 1  ✖ ℹ 🗑

**Scale** ℹ

Custom Alt-R™ gRNA, 2 nmol ▾

**Sequence** ✖ (5' → 3')

5' MOD ▾ INTERNAL ▾ 3' MOD ▾ BASES ▾

/AltR1/rArU rGrCrA rUrArG rCrUrA rGrArC rUrArG  
rArUrC rArGrA rCrUrC rUrCrU rCrGrA rUrCrA

# Bases: 38 (Min:30 Max:150)

GC: 44.7% Tm: 64.8°C ⚙ DeltaG: -62.25 kcal/mole

CONVERT TO RNA

Step 1: Enter Sequences (1 item)

Step 2: Order Custom CRISPR Essentials (0 items)

CONTINUE >

Show CRISPR Help

- Online ordering of any custom gRNA: including Cas13, pegRNA, gRNAs for novel nuclease systems
- Lengths: 30 to 150 nt
- Scales: 2, 10, 50, & 100 nmols
- Supports: RNA bases, Alt-R modifications, 2'OMe bases, and phosphorothioate linkages

# IDT CRISPR PROTEIN ENGINEERING

- Cas9 vs. Cas12a (Cpf1)
- Alt-R WT Cas9 vs. **Alt-R HiFi Cas9**
  - Problem with Cas9: off-target editing
  - Bacterial mutagenesis to evolve high-fidelity Cas9
- Alt-R WT Cas12a (Cpf1) vs. **Alt-R Cas12a (Cpf1) Ultra**
  - Problem with Cpf1: low on-target editing
  - Bacterial mutagenesis to evolve enhanced activity A.s. Cas12a
  - Actively working on developing mutants of *L.b.* Cas12a with increased activity
- Fluorescent CRISPR fusion proteins:
  - **Alt-R™ S.p. Cas9-GFP V3, 100 µg**
  - **Alt-R™ S.p. Cas9-RFP V3, 100 µg**

IDT-exclusive mutants

Confidential - Company Proprietary


# NUCLEASE

Product Name	Catalog #
Alt-R™ S.p. Cas9 Nuclease V3, 5 mg	10000735
Alt-R™ A.s. Cas12a (Cpf1) Ultra, 5 mg	10007804
Alt-R™ S.p. HiFi Cas9 Nuclease V3, 5 mg	10007803
Alt-R™ Cas9 Elec Enhancer, 100 nmol	10007805
Alt-R™ S.p.Cas9 V3, glycerol-free, 100µg	10007806
Alt-R™ S.p.Cas9 V3, glycerol-free, 500µg	10007807
Alt-R™ S.p.Cas9 V3, glycerol-free, 5 mg	10007808
Alt-R™ L.b. Cas12a (Cpf1) Ultra, 100 µg	10007922
Alt-R™ L.b. Cas12a (Cpf1) Ultra, 500 µg	10007923
Alt-R™ L.b. Cas12a (Cpf1) Ultra, 5 mg	10007924

## Large Scale Cas9 Proteins:


- Glycerol-free WTCas9v3 – ideal for sensitive cell lines
- As and Lbcas12a (Ultra) – useful for targeting AT-rich regions without available Cas9-specific PAM sequences. Much higher on-target potency than wild-type A.s. Cas12a (Cpf1). Can recognize many TTTT PAM sites in addition to TTTV motifs, increasing target range. Active at room temperature, making them flexible tools for applications requiring delivery at lower temperatures.

# A.s. Cas12a *ULTRA* HAS EDITING ACTIVITY COMPARABLE TO *S.p.* Cas9



71 Cas9/Cas12a RNP (TTTV PAMs) targeting multiple loci delivered by electroporation in HEK293 and Jurkat cells

# A.s. AND L.b. Cas12a *ULTRA* HAVE INCREASED TEMPERATURE TOLERANCE

# IDT FLUORESCENT CRISPR PROTEINS RETAIN HIGH ON-TARGET ACTIVITY



- *S.p.* Cas9 (Cas9 V3, Cas9-GFP, or Cas9-mCherry) was delivered at 2.0  $\mu$ M RNP into HEK293 cells by Nucleofector™ Nucleofection (Lonza) targeting sites within the HPRT gene
- *A.s.* Cas12a RNP (Cas12a *Ultra* or Cas12a-GFP) was delivered at a suboptimal dose of 50 nM RNP into HEK293 cells by nucleofection to achieve a range of editing across the tested HPRT sites



# Cas9-GFP ALLOWS FOR VISUALIZATION OF DELIVERY

2  $\mu$ M RNP delivery by Nucleofector™ Nucleofection (Lonza) into HEK293 cells, imaged ~18 hours after delivery using a confocal microscope



**Cas9  
V3**

**Cas9-  
GFP**


# DELIVERY PROFILE IS REFLECTIVE OF EDITING


# FLUORESCENTLY-LABELED tracrRNA

- Labeled tracrRNA was duplexed with Alt-R™ crRNA XTs
- crRNA:tracrRNA duplex complexed with Alt-R™ Cas9 V3
- RNP was delivered into cells using Lipofectamine™ RNAiMAX (Thermo Fisher) at 10 nM
- HEK293 cells were incubated for ~18 h at 37°C.
- Cells were washed twice with PBS and run on a flow cytometer


Cas9 V3  
+ Unlabeled  
tracrRNA

Cas9 V3  
+ Labeled  
tracrRNA

# ENRICHMENT OF EDITED CELLS USING LABELED tracrRNA BY FACS

High: Top 20%, Medium 80–60%, Low: Bottom 60% of cells based on fluorescent signal

(10 nM RNP delivery by Lipofectamine™ RNAiMAX (Thermo Fisher) into HEK293 cells, sorted after ~18 hours)

# IDT SOLUTIONS TO IMPROVE HOMOLOGY-DIRECTED REPAIR

---

# USE THE Alt-R HDR DESIGN TOOL FOR OPTIMAL DESIGN

## Alt-R HDR Design Tool

Design and order homology-directed repair (HDR) donor templates and associated Cas9 guide RNAs for genome editing human, mouse, rat, zebrafish, or C. elegans targets.

[← BACK](#) **1. PICK TARGET(S)** **2. SPECIFY DESIGN** [WATCH VIDEO ▶](#)

Visualize & specify your mutation by opening the sliders on the interactive map or input fields below. Need additional assistance? Click 'Watch video' to learn more.

NCBI Transcript Accession: NM\_000518.4

Exon Intron guide RNA Reference Sequence SNP/MNP Insertion Deletion

5226994 5227026

5' 3'

Start edit: 5226994 SNP/MNP/Deletion length: 0 [ZOOM TO EDIT](#) [PREVIEW TRANSLATION](#) [← PREVIOUS](#) [NEXT →](#)

Mutation:  [Select tag \(optional\)](#) [INSERT](#)

Homology arms: Left  Right   Use default homology

Number of designs: 5 Add silent mutations:  No  Yes

Guide: GTAACGGCAGACTTCTCCTC

[← BACK](#) [DESIGN >](#)


For S.p. Cas9 designs, guidelines based on empirical data are used to recommend donors that are associated with guide RNAs that are near the desired mutation and have both high on-target and off-target scores. For S.p. Cas9 D10A nickase designs, guidelines based on empirical data were used to recommend paired guide RNAs that flank the desired mutation, have PAM sites oriented away from the mutation, and have cut sites 37–70 nucleotides apart.

## Alt-R HDR Design tool: a novel bioinformatics tool for ssDNA HDR template design

- Human, mouse, rat, zebrafish, nematode, or custom input
- gRNA selection using IDT Alt-R gRNA design tool
- Addition of silent mutations improves rates of HDR
- Supports WT and Cas9 Nickase (D10A) strategies
- Single or multiple (batch) design
- Alt-R modified donors available in output

# HDR ENHANCER PROTEIN

## Alt-R HDR Enhancer V2

CRISPR-induced double-strand breaks are repaired via NHEJ or HDR. The Alt-R™ HDR Enhancer Protein promotes HDR by inhibiting 53BP1, making it specific to the HDR pathway, in contrast to DNA-PK inhibitors that enhance HDR by blocking the NHEJ pathway directly.

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
- ① ssDNA donor with ends homologous to DSB
- ② dsDNA donor with ends homologous to DSB
- ③ Plasmid donor with ends homologous to DSB

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


Can modifications reduce the risk of blunt integration with dsDNA?

# SOLUTIONS TO REDUCE BLUNT INTEGRATION WITH dsDNA HDR TEMPLATES


# Alt-R HDR Donor BlocksとAlt-R HDR Enhancer V2の併用が large HDR insertions効率を改善


*Lonza Nucleofector™ delivery system:*

K562 or HEK-293 cells, 2  $\mu$ M Cas9 RNP, 50 nM ss or dsDNA (GFP tag, 200bp HA), 0 or 1  $\mu$ M Alt-R HDR Enhancer V2 (24 hrs) (n=4)

# HDR ENHANCER V2 ALSO IMPROVES HDR RATES FOR LARGE INSERTIONS

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
## rhAmpSeq™ MULTIPLEXED AMPLICON SEQUENCING FOR CRISPR

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# EXAMPLE CRISPR EXPERIMENT WORKFLOW

# rhAmp PCR TECHNOLOGY INNOVATED AT IDT

*BMC Biotechnology* (2011) 11:80

METHODOLOGY ARTICLE

Open Access

## RNase H-dependent PCR (rhPCR): improved specificity and single nucleotide polymorphism detection using blocked cleavable primers


Joseph R Dobosy, Scott D Rose, Kristin R Beltz, Susan M Rupp, Kristy M Powers, Mark A Behlke<sup>\*</sup> and Joseph A Walder

Blocked inactive primers anneal to target

RNase H2 recognizes hybridized internal RNA base

RNase H2 cleaves hybridized primers

DNA polymerase extends newly unblocked primers

# rhAmpSeq DESIGN TOOL

SARS-CoV-2. Don't let up. We'll help.

GET HELP

PRODUCTS & SERVICES ▾

APPLICATIONS & TECHNOLOGIES ▾

SUPPORT & EDUCATION ▾

TOOLS ▾

COMPANY ▾

ORDER MENU

You are using the IDT portal.

## rhAmpSeq Design Tool

Welcome to the IDT rhAmpSeq design tool. If you encounter any questions while using the tool, contact [applicationsupport@idtdna.com](mailto:applicationsupport@idtdna.com). If you don't see what you're looking for, submit a request to the rhAmpSeq Custom Design Service here.

Create new design

Review design results

Dashboard

Application

Input format

Species

Minimum insert size

Maximum insert size

Panel name

DESIGN

CLEAR AND RESET ↶

[Open tool resources](#)

[Open user guide](#)

Paste/Type input

Upload file

Enter your coordinates separated by tab below.

chr7 12345 12365

Proprietary


# rhAmpSeq DESIGN TOOL

## Supported Reference Genomes

Species	Reference genome
<i>Homo sapiens</i> (human)	GRCh38
<i>Homo sapiens</i> (human)	GRCh37
<i>Mus musculus</i> (mouse)	GRCm38
<i>Rattus norvegicus</i> (rat)	Rnor6.0
<i>Bos taurus</i> (cow)	ARS UCD1.2
<i>Bos taurus</i> (cow)	UMD3.1
<i>Brassica napus</i> (canola)	Bra napus v2.0
<i>Caenorhabditis elegans</i> (nematode)	WBcel235
<i>Danio rerio</i> (zebrafish)	GRCz11
<i>Gallus gallus</i> (chicken)	GRCg6a
<i>Glycine max</i> (soybean)	Glycine max v2.1
<i>Gossypium hirsutum</i> (cotton)	ASM98774v1
<i>Oryza sativa</i> (rice)	IRGSP 1.0
<i>Solanum lycopersicum</i> (tomato)	SL3.0
<i>Solanum tuberosum</i> (potato)	SolTub 3.0
<i>Sus scrofa</i> (pig)	Sscrofa11.1
<i>Triticum aestivum</i> (wheat)	IWGSC v1.0
<i>Zea mays</i> (maize)	B73 RefGen v4

# rhAmpSeq TECHNOLOGY AND Alt-R HiFi Cas9

Verification and validation Alt-R *S.p.* HiFi Cas9 Nuclease V3

## rhAmpSeq system savings:

- 8 x 40 assay designs
- ~1500 individual PCRs reduced to <96
  - Master mix
  - gDNA
- Library quantification
- Full-time equivalent hours
  - Months to days

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CRISPR with confidence. 詳しくはこちらから。

### Aldevron SpCas9 Nuclease and SpyFi® Nuclease

ADD TO ORDER

Recombinant *Streptococcus pyogenes* Cas9 nuclease, purified from an *E. coli* strain expressing the SpCas9 nuclease. Contains N- and C-terminal SV40 nuclear localization sequences (sNLS) and no affinity purification tags. Provided in solution at 10 mg/mL. SpyFi® Cas9 nuclease has similar on-target potency to wild-type SpCas9, but with reduced off-target effects, allowing for more precise genome editing.

Quantity	Product	Catalog #	Price
<input type="text" value="0"/>	sNLS-SpCas9-sNLS Nuclease, 0.25 MG	10017687	¥44,760 JPY
<input type="text" value="0"/>	sNLS-SpCas9-sNLS Nuclease, 5 MG	10017688	¥744,600 JPY
<input type="text" value="0"/>	SpyFi Cas9 Nuclease, 0.25 MG	10017689	¥69,700 JPY
<input type="text" value="0"/>	SpyFi Cas9 Nuclease, 5 MG	10017690	¥946,500 JPY


### Aldevron Eureka-V™ Nuclease

ADD TO ORDER

Recombinant Class 2, Type V CRISPR-Cas nuclease, purified from an *E. coli* strain expressing the recombinant nuclease. Contains N- and C-terminal SV40 nuclear localization sequences (sNLS) and no affinity purification tags. Provided in solution at 10 mg/mL. Requires acknowledgment of Limited Use Label License (LULL).


Quantity	Product	Catalog #	Price
<input type="text" value="0"/>	Eureka-V Nuclease, 1MG	10017800	¥211,500 JPY
<input type="text" value="0"/>	Eureka-V Nuclease, 5MG	10017811	¥964,700 JPY

# High-efficiency multiplexed T cell editing

**High multiplexed knockout efficiency with wild type (WT) and SpyFi RNP.** To knockout both *TCR* and *PD-1*, activated human T cells were electroporated with an equimolar mixture of RNP targeting *TRAC*, *TRBC*, and *PD-1*. Expression was measured by flow cytometry 5 days post electroporation. Expression was substantially reduced following electroporation with WT or SpyFi RNPs compared to no electroporated cells,  $n = 1$ .

## Reduced off-target effects

**Figure 2. Fewer off-target events with *SpyFi* vs WT RNP.** SpCas9 or *SpyFi* Cas9 RNP targeting the *TRAC* locus was transfected into activated human T cells using Maxcyte electroporation (EP) technology. The frequency of off-targeting editing for each RNP was determined by measuring InDel percentages at four off-target binding sites previously identified by iGUIDE [1]. DNA flanking off-target binding sites in *CLIC2*, *LDAH*, *ANKS1B*, and *ADCY10* were amplified by PCR, gel-purified, and sequenced by next generation sequencing (NGS). (A) At the *CLIC2* locus *SpyFi* RNP reduced off target editing by 94% compared to WT RNP. (B) Similarly, at the *LDAH* locus there was an 87% reduction in off-target editing. At both these loci *SpyFi* had off-target effects comparable to no electroporation controls. At two other sites (*ANKS1B* and *ADCY10*), off-target editing was generally less frequent but was diminished by *SpyFi* (33% less at *ANKS1B* locus, 15% less at *ADCY10*, data not shown),  $n=1$ .

# CRISPR cGMP gRNA Manufacturing Overview

For customers looking to accelerate their CRISPR gene editing program from discovery to clinical trials, our Engineering Run and full cGMP Compliant gRNA Manufacturing Services offer a streamlined and regulated solution. Backed by comprehensive documentation, our guide RNA CRISPR services simplify regulatory filings, offering a straightforward path to clinical success.

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	Engineering Run	cGMP
Availability	Delivery times after placing the order	Delivery times after placing the order
Cleanroom	Certification not required	ISO 8 Clean Room - Certified
Changeover	In-place with cleaning validation	In-place with cleaning validation
Materials	QA release on raw materials	QA release on raw materials
Batch Records	Draft based on customer's specifications	Based on customer specification
Release Testing	Qualified methods	Validated methods

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